

Arial® Pendant & Gateway with Auto Fall Detect

OVERVIEW



Receive an automatic alert when a resident may have fallen



ABOUT

Arial Pendants offer a comfortable and discreet way for residents to call for help anywhere in your community. New, advanced Auto Fall Detect technology is designed to automatically send an alert when a fallen resident is unable to press the pendant call button for help. The Auto Fall Detect feature is designed to help detect a fall and send an alert to the Arial system, notifying caregivers to respond to a potential fall event.

HOW IT WILL WORK

Resident falls > Algorithm confirms > A unique fall alarm is sent to Arial system via cloud



Device sensors look for a height change >2.5 ft., rotation of >30° and acceleration

Thirty second calculation window to rule out a recovery and/or sensor fluctuations

The pendant vibrates, its red light flashes, and staff are alerted to respond to a possible fall event

4600 Vine Street
Lincoln, NE 68503
+1 888 622 6992

North America
info@securitashealthcare.com

International
global@securitashealthcare.com

securitashealthcare.com

About Securitas Healthcare

Securitas Healthcare empowers caregivers to deliver connected, productive and safe care. Its innovative portfolio of healthcare solutions helps over 15,000 hospitals, clinics and senior living organizations worldwide protect people, use assets efficiently and understand their operations for a caring and healing environment. Securitas Healthcare is proud to be part of Securitas, the world's leading intelligent protective services partner. For more information, visit us at securitashealthcare.com

DOC-12-50042-AB

Arial is a trademark of Securitas Healthcare LLC and/or its affiliates. © 2023 Securitas Healthcare LLC. All rights reserved. This publication is accurate as of the date of printing.